



### Key Features

- Refrigerator & Freezer Combination
- Strong Refrigeration
- Five Alarm Functions

### More Detail:

[http://www.haiermedical.com/products/preservation/combined\\_refrigerator\\_and\\_freezer/hycd-282.html](http://www.haiermedical.com/products/preservation/combined_refrigerator_and_freezer/hycd-282.html)

## Specifications

### Technical Data

Type	Combined Refrigerator and Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Refrigerator: Forced air cooling Freezer: Direct cooling
Defrost Mode	Manual/Auto
Refrigerant	CFC-Free
Noise Level (dB)	43

### Performance

Temperature Range( )	2~8/-20~-40
----------------------	-------------

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	400

<b>Electrical Current(A)</b>	2.8
<b>Dimensions</b>	
<b>Capacity(L/Cu.Ft)</b>	Refrigerator: 185/6.5 Freezer: 97/3.4
<b>Net/Gross Weight(approx)</b>	145/160(kg) 320.0/352.7(lbs)
<b>Interior Dimensions(W*D*H)</b>	Refrigerator: 605*510*720(mm) 23.8*20.0*28.3(in) Freezer: 515*465*440(mm) 20.3*18.3*17.3(in)
<b>Exterior Dimensions(W*D*H)</b>	736*660*1810(mm) 29.0*26.0*71.3(in)
<b>Packing Dimensions(W*D*H)</b>	875*770*1920(mm) 34.4*30.3*75.6(in)
<b>Container load (20'/40'/40'H)</b>	12/26/26
<b>Functions</b>	
<b>High/Low Temperature</b>	Y
<b>Remote Alarm</b>	Y
<b>Power Failure</b>	Y
<b>Sensor Error</b>	Y
<b>Low Battery</b>	Y
<b>High Ambient Temp</b>	Y
<b>Accessories</b>	
<b>Caster</b>	Y
<b>Foot</b>	Y
<b>Test Hole</b>	Y
<b>Shelves/Inner Doors</b>	3/-
<b>Drawers</b>	2
<b>USB Interface</b>	Y
<b>Others</b>	
<b>Certificate</b>	CE