

Liquid Nitrogen Container For Storage

Products

Features

Model	Unit	YDS-20B	YDS-30B	YDS-30B-80	YDS-35B	YDS-35B-80	YDS-35B-125	
Capacity	L	20	31.5	31.5	35	35	35	
Empty Weight	Kg	11.5	12.5	13	14.1	14.3	14.8	
Neck Diameter	mm	50	50	80	50	80	125	
Outer Diameter	mm	394	461	461	461	461	461	
Height	mm	656	656	656	695	698	700	
Static Evaporation Loss	L/d	0.20	0.19	0.28	0.20	0.29	0.41	
Static Hold Time	d	101	159	106	179	119	86	
Canister Outer Diameter	mm	38	38	63	38	63	97	
Canister Height	mm	120/276	120/276	120/276	120/276	120/276	120/276	
Numbers of Canister	EA	6(3)	6(3)	6(3)	6(3)	6(3)	6(3)	
Straw Holding Capacity(single level)	0.5ml	EA	792	792	2244	792	2244	5124
	0.25ml	EA	1788	1788	5022	1788	5022	11952
Straw Holding Capacity(single level)	0.5ml	EA	1284	1284	3624	1284	3624	9048
	0.25ml	EA	2832	2832	8460	2832	8460	19944

Model	Unit	YDS-50B	YDS-50B-80	YDS-50B-125			
Capacity	L	50	50	50			
Empty Weight	Kg	19.3	20.1	21			
Neck Diameter	mm	50	80	125			
Outer Diameter	mm	461	461	461			
Height	mm	770	766	753			
Static Evaporation Loss	L/d	0.23	0.34	0.45			
Static Hold Time	d	213	147	110			
Canister Outer Diameter	mm	38	63	97			
Canister Height	mm	120/276	120/276	120/276			
Numbers of Canister	EA	6(3)	6(3)	6(3)			
Straw Holding Capacity(single level)	0.5ml	EA	792	2244	5124		
	0.25ml	EA	1788	5022	11952		
Straw Holding Capacity(single level)	0.5ml	EA	1284	3624	9048		
	0.25ml	EA	2832	8460	19944		